

Journalbeiträge

1. Chapochnikov NM, Takago H, Huang CH, Pangršič T, Khimich D, Neef J, Auge E, Göttfert F, Hell SW, Wichmann C, Wolf F, Moser T (2014) Uniquantal release through a dynamic fusion pore is a candidate mechanism of hair cell exocytosis. *NEURON* 83(6): 1389-403, doi: 10.1016/j.neuron.2014.08.003
2. Goßler C, Bierbrauer C, Moser R, Kunzer M, Holz K, Pletschen W, Köhler K, Wagner J, Schwaerzle M, Ruther P, Paul O, Neef J, Keppeler D, Hoch G, Moser T, Schwarz U (2014) GaN-based micro-LED arrays on flexible substrates for optical cochlear implants. *J PHYS D APPL PHYS* 47: 205401-205406, doi: 10.1088/0022-3727/47/20/205401
3. Hernandez VH, Gehrt A, Jing Z, Hoch G, Jeschke M, Strenzke N, Moser T (2014) Optogenetic stimulation of the auditory nerve. *J Vis Exp* 1(92): e52069, doi: 10.3791/52069
4. Hernandez VH, Gehrt A, Reuter K, Jing Z, Jeschke M, Mendoza Schulz A, Hoch G, Bartels M, Vogt G, Garnham CW, Yawo H, Fukazawa Y, Augustine GJ, Bamberg E, Kügler S, Salditt T, de Hoz L, Strenzke N, Moser T (2014) Optogenetic stimulation of the auditory pathway. *J CLIN INVEST* 124(3): 1114-29, doi: 10.1172/JCI69050
5. Mendoza Schulz A, Jing Z, Sánchez Caro JM, Wetzel F, Dresbach T, Strenzke N, Wichmann C, Moser T (2014) Bassoon-disruption slows vesicle replenishment and induces homeostatic plasticity at a CNS synapse. *EMBO J* 33(5): 512-27, doi: 10.1002/embj.201385887
6. Neef J, Jung S, Wong AB, Reuter K, Pangrsic T, Chakrabarti R, Kügler S, Lenz C, Nouvian R, Boumil RM, Frankel WN, Wichmann C, Moser T (2014) Modes and regulation of endocytic membrane retrieval in mouse auditory hair cells. *J NEUROSCI* 34(3): 705-16, doi: 10.1523/JNEUROSCI.3313-13.2014
7. Revelo NH, Kamin D, Truckenbrodt S, Wong AB, Reuter-Jessen K, Reisinger E, Moser T, Rizzoli SO (2014) A new probe for super-resolution imaging of membranes elucidates trafficking pathways. *J CELL BIOL* 205(4): 591-606, doi: 10.1083/jcb.201402066
8. Weiler S, Krinner S, Wong AB, Moser T, Pangršič T (2014) ATP hydrolysis is critically required for function of CaV1.3 channels in cochlear inner hair cells via fueling Ca²⁺ clearance. *J NEUROSCI* 34(20): 6843-8, doi: 10.1523/JNEUROSCI.4990-13.2014
9. Wong AB, Rutherford MA, Gabrielaitis M, Pangrsic T, Göttfert F, Frank T, Michanski S, Hell S, Wolf F, Wichmann C, Moser T (2014) Developmental refinement of hair cell synapses tightens the coupling of Ca²⁺ influx to exocytosis. *EMBO J* 33(3): 247-64, doi: 10.1002/embj.201387110